

What is Buzzy®?

Buzzy® is an award-winning FDA-cleared medical device to CONTROL PAIN.



-  MOST PROVEN
-  FAST
-  EASY
-  REUSABLE
-  COST EFFECTIVE

How Does Buzzy® Work?

Using a combination of vibration and cold to desensitize nerves, Buzzy® controls pain on contact.




Reduce Pain Up To 88%



PROVEN in 60+ clinical trials.
IV access 45 minutes faster.
Saves ~\$13.00 per point of pain reduction.

Using Buzzy® Is Easy



Secure frozen wings on vibration unit.

"Buzz the skin to sleep" for one minute.

Do procedures with Buzzy® between brain & pain.



Search & find cards clinically proven to reduce pain!

- Snaps onto a name badge or key ring.
- Age-appropriate games & tasks on the back.
- One infant card per set.

Accessories available for sale on Shop.PainCareLabs.com



"I had no idea that she had even stuck me...This is a huge step in life for me."

- Matt S.

"My father was going to quit dialysis until we found Buzzy®. Dr. Baxter, you're a lifesaver!"

- Karen S.

Model Comparison



Personal

- Bug bites & bumps, home injections, diabetes management, lab draws.
- Toggle switch for continuous vibration.
- Includes 2 reusable single-user reversible ice wings, 2 AAA batteries, and hands-free strap.



Healthcare XL

- Longer procedures, infusions, IVs.
- Internal tourniquet slot.
- Toggle switch for continuous vibration.
- Includes 4 reusable multi-patient ice wings, 2 AAA batteries.



Healthcare Mini

- Vaccinations, small hands, heel & finger sticks.
- Automatic 3-minute shutoff.
- Includes 4 reusable multi-patient ice wings, 2 AAA batteries.

NEW Buzzy® Pro



Dialysis • Adult Lab Procedures

PainCareLabs.com



POWER OVER PAIN



Home Procedures
IVF, growth hormone, insulin,
methotrexate, Humira®,
glucose monitoring, Lovenox®



Dermatology
Lidocaine & steroid injections,
Dupixent®, Botox®, *C. albicans*,
hair removal, electrolysis



Clinic Use
IV access, vaccinations,
dialysis, lidocaine injection,
allergy testing, blood draws



Post-surgery, ACL injury,
wrist pain, plantar fasciitis
knee pain, oral surgery,
carpal tunnel, & more

Coming Soon - Low Back Pain Relief



Invented by a physician, Buzzy® is FDA cleared to control pain associated with injections, venipuncture, IV starts, cosmetic injections, and the temporary relief of minor injuries (muscle or tendon aches, splinters, insect stings).

Pain Care Labs
195 Arizona Ave NE, LW08
Atlanta, GA 30307 USA
877-805-BUZZ(2899)
info@PainCareLabs.com



The Best Solution for Needle Pain



"Buzzy has literally changed my life!
Needle phobia GONE!"
- Angela Vermillion

PainCareLabs.com

Rev. 09/23/2021
FDAClearance
Registered (510k)
CLEARED

